

QuadraPure™ SP Cartridges

The QuadraPure™ range of metal scavengers offers practical methods for removal of metal catalyst residues, allowing easier, faster and cleaner processes to be developed. In addition to the loose resin, QuadraPure™ TU is also conveniently available in cartridge format compatible with Biotage SP automated chromatography units.



QuadraPure™ TU SP Cartridges

QuadraPure™ Cartridge Advantages:

| | |
|----------------------------------|---|
| Cleaner products | typically achieves less than <1-10 ppm residual metal contamination |
| Fast, efficient processes | selective metal removal in single pass through cartridge |
| Ease of handling | avoids handling/filtering loose beads |
| High selectivity | minimal product absorption |
| Process intensification | higher metal loadings with faster kinetics than in batch |
| Improved yields | versus conventional work-up and isolation steps |
| Process development | scalable to manufacturing scale cartridges |

QuadraPure™ SP Cartridge Characteristics:

| Cartridge size | Approx mass of QuadraPure™ | Approx cartridge metal capacity | Optimum flow rate* |
|----------------|----------------------------|---------------------------------|--------------------|
| 12M | 8 g | 1 mmol | 1 ml/min |
| 25S | 18 g | 2 mmol | 2 ml/min |
| 40M | 80 g | 10 mmol | 8 ml/min |

*Optimum flow rate is ~5 column volumes/h


For pricing and ordering of QuadraPure™ cartridges please contact info@reaxa.com




QuadraPure™ SP Cartridges

Examples of QuadraPure™ TU SP Cartridge Scavenging Applications:

1. Efficient removal of multi-valent mixed metals in the presence of ligand and base following a Sonogashira reaction:




Cu 54 ppm, Pd 58 ppm
PPh₃ & Et₃N in crude solution



Cu and Pd <1 ppm,
zero product absorption

2. Highly selective removal of palladium species, compared to standard laboratory methods, following a Suzuki reaction and work up:



| Method | % Pd removed | % product absorbed |
|--------------------------------|--------------|--------------------|
| Silica (pad) | 6 | 2 |
| Carbon cartridge | 56 | 28 |
| QuadraPure™ TU cartridge (12M) | >99 | <1 |

| Product | Structure | Metals removed | Availability |
|------------------|-----------|---|---|
| QuadraPure™ TU | | Pd, Pt, Ru, Rh, Au, Ag, Cu, Hg, Pb, Cd, Ni, Co, Fe, V, Zn | 12M, 25S ex stock 40M info@reaxa.com |
| QuadraPure™ IDA | | Cu, Al, Ga, In, V, Pb, Ni, Zn, Cd, Be, Mn, Sr, Ba, Co, Fe | info@reaxa.com |
| QuadraPure™ AMPA | | Fe, Cu, Ni, Al, Co, V | info@reaxa.com |
| QuadraPure™ BZA | | Rh, Pd, Cu, Co, Ni and electrophiles* | info@reaxa.com |

* Literature example: A Fully Automated Continuous Flow Synthesis of 4,5-Disubstituted Oxazoles, M. Baumann, IR. Baxendale, S. V. Ley, C. D. Smith, G. K. Tranmer, *Org Lett.*, 2006, 8, 5231.

For pricing and ordering of QuadraPure™ cartridges please contact info@reaxa.com

